

New!
AERO/VENT™
"Jr"

Radioaerosol Lung Imaging System
Model AV-100 Series

The power of two.

- **Two tube system**
(increased speed with *no loss* of image quality)
- **Two convenient models**
(pole-mount and table top)
- **Two times as fast***
(1.5-2.5 minute full dosing)
- **Two filters available**
(Bacteria and HEPA)

No charge loaner shield available.

* as the Aero/Vent™ Plus



AV-101 Pole Mount Shield



AV-102 Table Top Shield

It's about quality.
It's about time.

Aero/Vent™ "Jr": Our most economical system!

AERO/VENT™ "Jr"

Features

- Portability and Safety
- Reduced Breathing Time
- Low breathing resistance
- Low deadspace
- Excellent images
- Useful accessories
- Expert customer support



Optional UV-203
Universal Cart



For over a decade and a half, Medi/Nuclear has been dedicated to producing the highest quality radioaerosol lung scanning systems in the field. All of its fine Insta/Vent™ and Aero/Vent™ products offer unmatched image quality as a result of a particle size of 0.3 µm (MMAD). With the same attention to safety, patient comfort and technologist convenience, the Aero/Vent™ "Jr" now offers technologists a very economical 2-tube breathing system that can be used with children as well as adults to reduce deadspace, increase patient comfort and reduce breathing time.

Ordering Information

Product Description	Cat. No.	Packaging
Aerosol Breathing Circuits :		
Aero/Vent™ "Jr" with Bacteria Filter & Rigid Mouthpiece	AV-100B	20/case
Aero/Vent™ "Jr" with Bacteria Filter & Mask	AV-100BM	20/case
Aero/Vent™ "Jr" with HEPA Filter & Rigid Mouthpiece	AV-100H	20/case
Aero/Vent™ "Jr" with HEPA Filter & Mask	AV-100HM	20/case
Shields and Accessories:		
Pole Mount Shield	AV-101	each
Table Top Shield	AV-102	each
Universal Cart	UV-203	each