

New!

AERO/VENTTM
MAX

Radioaerosol Lung Imaging System
Model AV-500 Series

It's about time.



- **Two tube system**
Increased speed with *no* loss of image quality
- **No-charge loaner shield available**
- **HEPA filter standard.**
- **Mask configuration option**
- **1.5-2.5 minute full dosing**

It's about quality.

Model AV-501 Shield

AERO/VENTTM MAX

Features

- Portability and Safety
- Reduced Breathing Time
- Low breathing resistance
- Low deadspace
- Excellent images
- Useful accessories
- Expert customer support



Optional UV-203
Universal Cart
Shown with Aero/Vent Jr



For over a decade and a half, Medi/Nuclear has been dedicated to producing the highest quality radioaerosol lung scanning systems in the field. All of its fine Insta/VentTM and Aero/VentTM products offer unmatched image quality as a result of a particle size of 0.3 μm (MMAD). With the same attention to safety, patient comfort and technologist convenience, the Aero/VentTM "Max" now offers technologists a very economical 2-tube breathing system that can be used with children as well as adults to reduce deadspace, increase patient comfort and reduce breathing time.

Ordering Information

Product Description	Cat. No.	Packaging
Aerosol Breathing Circuits:		
Aero/Vent TM "Max" with HEPA Filter & Safety Shield Mouthpiece	AV-500 B	20/case
Aero/Vent TM "Max" with HEPA Filter & Mask	AV-500 M	20/case
Shields and Accessories:		
Table Top Shield	AV-501	each
Universal Cart	UV-203	each